



FLEXXIDISC

Operating Instructions

HEADQUARTERS
Horse Dental Equipment
Châteaubourg (35)
FRANCE

Ph: +33 (0)2 99 00 71 29
www.horse-dental-equipment.com

AMERICA
HDE North America Ltd.
Houston (TX)
USA

Ph: (713) 444-2531
www.horse-dental-equipment.us

www.horse-dental-equipment.com

Precautions

- Maintenance must be done at least once a year or every 500 horses
- Do not use this handpiece for something else than its intended purpose
- Disassembly of the flexxidisc's head will void the warranty

The different models

MODEL	CURVED TWIST	STRAIGHT	LENGTH (CM)
LONG	✓	✓	54
STANDARD	✓ <small>(with or without extension)</small>	✓ <small>(with or without extension)</small>	32
SHORT		✓	18


FOR ANY PROBLEM, PLEASE CONTACT
THE OFFICIAL TECHNICAL CENTER

DIFFERENT SHAPES & SIZES
ADAPTED TO ALL DENTAL PRACTICES



www.horse-dental-equipment.com

Step n°1



Hold both head + ring and pull **strongly** outward.

Step n°2



Turn gently and release.

Step n°3



The head is correctly positioned when it reaches the notch and you hear a click.

Diamond burrs

1. Start the motor with the handpiece attached
2. Run it at high speed with the head pointing down to remove any remaining water
3. Clean the diamond burr by running it over the HDE stainless steel brush
4. Continue until the remains of dentine have completely disappeared
5. Clean with a cloth moistened with clean water and a disinfectant
6. Thoroughly dry with a clean, soft cloth
7. Store the burr in its dedicated place within your clean, dry case

The burr is resistant to water projection, do not fully immerse in water
Always use clean water without corrosive nor abrasive products
Only use the HDE stainless steel brush



Drive shafts

1. Clean the outside sheath with a damp sponge
2. Dry thoroughly with a clean, soft cloth
3. When finished working, detach the drive shaft from the motor and the handpiece
4. Store the drive shaft in your clean, dry case

Ensure not to bend the drive shaft in an excessive radius of curvature

Motor

1. Turn off the motor before cleaning
2. Clean the motor with a damp sponge
3. Thoroughly dry the motor with a clean, soft cloth
4. Store the motor in your clean, dry case

Case

1. Clean the case with a damp sponge
2. Thoroughly dry the case with a clean, soft cloth
3. Store in a clean and well ventilated place

Make sure the case is opened from time to time, to ensure your equipment is well ventilated

Advantage n° 1



An efficient floating of the sharp enamel points thanks to its diamond disc.

Advantage n° 2



The very slim head & body profile allows an easy access, far into the mouth.

The curved shape allows to see the disc working on the teeth, especially when working on the vestibular part.

Advantage n° 3



Its rubber handle offers optimized comfort and grip.