**HORSE DENTAL
EQUIPMENT**

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
MP 144350

Battery assembly

In order to operate properly, individual Li-Ion cells are mechanically and electrically integrated in battery assemblies specific to each application. The battery assembly incorporates electronics for performance, thermal and safety management.

Technology

- Graphite-based anode
- Lithium Cobalt oxide-based cathode
- Electrolyte: organic solvents
- Built-in redundant safety protections (shutdown separator, circuit breaker, safety vent)
- Batteries assembled from MP cells feature an electronic protection circuit

Built-in protection devices ensure safety in case of:

- Exposure to heat
- Exposure to direct sunlight for extended periods of time
- Short circuit
- Overcharge
- Overdischarge

When handling Saft MP batteries:

- Do not solder directly to cell terminal
- Do not disassemble
- Do not remove the protection circuit
- Do not incinerate

Transportation and storage:

- Store in a dry place at a temperature preferably not exceeding 30°C
- For long-term storage, keep the battery within a (30 ± 15) % state of charge

Saft
Specialty Battery Group

www.saftbatteries.com

Doc. N° 54044-2-1009




Information in this document is subject to change without notice and becomes contractual only after written confirmation by Saft.

Published by the Communications Department.

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Produced by Arthur Associates Limited.

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
Rechargeable lithium-ion battery MP 144350

High performance
Medium Prismatic cell

**Saft always supplies MP cells
in assemblies or in customized
battery system constructions**

Benefits

- Extended autonomy and life for mobile systems
- A broad operating temperature range
- Recommended for ruggedized designs
- Easy integration into compact and light systems
- Very high energy density (344 Wh/l and 143 Wh/kg)
- Unrivaled low temperature performance

Key features

- Excellent charge recovery after long storage, even at high temperature
- Maintenance-free
- Non-restricted for transport/ Non-assigned to Class 9 according to the UN Recommendations on the transport of dangerous goods - Model Regulations
- Underwriters Laboratories Component Recognition (File Number MH 12609) (UL)

Main applications

- Mobile asset tracking
- Small UPS
- Soldier of the future equipment
- Portable radios
- Professional portable lighting
- Bar code readers
- Portable payment terminal

Electrical characteristics

Nominal voltage (0.5 A rate at 20°C)	3.75 V
Typical capacity (at 0.5 A 20°C 2.5 V cut-off)	2.6 Ah
Nominal energy	9.75 Wh

Mechanical characteristics (Sleeved 100 % charged cell)

Thickness	14.9 mm
<i>(Thickness tends to increase with cycling, typically obtained after 500 cycles. Consult Saft) (At beginning of life 14.6 mm)</i>	
Width (max)	43.9 mm
Height (max)	54.5 mm
Typical weight	68 g
Lithium equivalent content	0.78 g
Volume	28 cm ³

Operating conditions

Charge method	Constant Current/Constant Voltage	
Charge voltage	4.20 +/- 0.05 V	
Maximum recommended charge current* *	2.6 A (C rate)	
Charge temperature range*	-20°C to +60°C	
Charge time at 20°C	To be set as a function of the charge current:	
	C rate	= 2 to 3 h
	C/2 rate	= 3 to 4 h
	C/5 rate	= 6 to 7 h
Maximum continuous discharge current* *	5.0 A (~2C rate)	
Pulse discharge current	up to 10 A (~4C rate)	
Discharge cut-off voltage	2.5 V	
Discharge temperature range*	-50°C to +60°C	

* For optimized charging below 0°C, 60°C and discharging at -50°C, consult Saft.

** Electronic protection circuits within battery packs may limit the maximum charge/discharge current allowable. Consult Saft.