

**HORSE DENTAL**  
**EQUIPMENT**

**Designed for Equine Dentistry**

[www.horse-dental-equipment.com](http://www.horse-dental-equipment.com)



HORSE DENTAL  
EQUIPMENT

# Equine Dentalscope





# Practitioner's benefits

- ✓ Clear and precise diagnosis
- ✓ Better patient follow-up
- ✓ Additional service and value
- ✓ Develops horse owners understanding



# Product characteristics



## General features

- ✓ 3 monitoring options:
  - TV screen
  - Wireless USB connector for phone / tablet / computer (optional)
  - Wireless screen (optional)
- ✓ 692 mm total length
- ✓ Head dimensions: width 17mm x thickness 17.2mm
- ✓ Rounded head on top
- ✓ Stainless steel
- ✓ Weight: 360 g

# Product characteristics



## Camera

- ✓ Video resolution: 704 x 576 pixels
- ✓ PAL video system
- ✓ 4 glass lenses
- ✓ 70° angle view
- ✓ 5.000 K led lighting

## Autonomy

- ✓ 1 hour (fully charged)



# Starter set



1 Dentalscope

2 Charger

3 1,5 m  
video cable

4 Storage case

# Starter set + wireless USB connector

## Monitor on Smartphone / Tablet / Computer

- Plug the radio frequency wireless connector on a USB socket, install the free DVR /OTG Software on the device (Smartphone, Tablet, Computer), and watch in live streaming.
- Dedicated software to take pictures and record videos
- Available on PC Win 7, Win 10, Android phones and tablets



# Starter set + wireless monitor

Display on a wireless radio frequency monitor



- 10" DVR wireless monitor with recording on SD-CARD (video + picture)
- 1 10" wireless monitor (1-hour autonomy fully charged)
- 2 Charger
- 3 Storage case







# Set options overview

	Starter set	Starter set + wireless USB connector	Starter set + wireless monitor
Dentalscope	✓	✓	✓
Charger	✓	✓	✓
Video cable (length: 1,5 m)	✓	✓	✓
Storage and transportation case	✓	✓	✓
USB connector + free software		✓	
10" DVR viewing monitor + micro-SD card + storage case			✓

# Special features



"Plug & Play"




Lightweight  
and portable




Sturdiness



Wireless  
options



Take clear  
pictures & videos



Battery inside  
the dentalscope

# Dedicated services

---



After-sales  
service




1 year  
warranty



# Video presentation

---



 Click on the picture to watch the full video

# Set prices

## Starter set



Dentalscope for connection to a display device with cable

From 2,499 €

## Starter set + wireless USB connector



Dentalscope with USB wireless connector for phone / tablet / computer

From 3,199 €


## Starter set + wireless monitor



Dentalscope with a 10" monitor including recording system, micro-SD card and storage case

From 3,799 €

\* EXW Châteaubourg (France), excluding transport costs, taxes and duties valid until 31 March 2022.



**HORSE DENTAL**  
**EQUIPMENT**

**Designed for Equine Dentistry**

[www.horse-dental-equipment.com](http://www.horse-dental-equipment.com)